# Rapid COVID-19 Test



- UK Government approved and tested
- COVID-19 antigen Rapid Test Device
- Results in <15 minutes

## Product Information

The Coronavirus Ag Rapid Test Cassette (Swab) is an in vitro immunochromatographic assay for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in direct nasopharyngeal (NP) swab specimens directly from individuals who are suspected of COVID-19 by their healthcare provider within the first ten days of symptom onset.

It is intended to aid in the rapid diagnosis of SARS-CoV-2 infections. Negative results from patients with symptom onset beyond ten days, should be treated as presumptive and confirmation with a molecular assay, if necessary, for patient management, may be performed.

The Coronavirus Ag Rapid Test Cassette (Swab) does not differentiate between SARS-CoV and SARS-CoV-2.

Information	In Detail
Time to result	15 minutes
Storage	2-30°C
Shelf life	24 months
Specimen type	Nasopharyngeal swab



### **Contents**

- 20 Test cassettes
- 20 Extraction tubes and tips
- 2 Extraction buffer vials
- 1 Workstation
- 1 Package insert

## Test Procedure & Interpretation

### **Specimen collection**


Use the nasopharyngeal swab provided in the kit


- 1. Carefully insert the swab in the patient's nostril.
- 2. Swab over the surface of the posterior nasopharynx and rotate the swab several times.
- 3. Withdraw the swab from the nasal cavity.




# Sample Preparation

1. Insert the test extraction tube into the workstation provided in the kit. Make sure that the tube is standing upright and reaches the bottom of the workstation.












4. Leave the swab in the extraction tube for 1 minute.




5. Squeeze the tube several times from the outside to immerse the swab. Remove the swab.




6. Push the tip which contains the filter onto the extraction tube. Ensure the tip has a tight fit.







# Test Procedure and Interpretation of Results



## Product Information

Clinical Performance of the Coronavirus Ag Rapid Test Cassette (Swab) was evaluated by being involved in 7 sites within the US, where patients were enrolled and tested. It has also passed phase 3 UK testing for government approval.

Testing was performed by 24 Health Care Workers that were not familiar with the testing procedure. A total of 317 fresh nasopharyngeal swab samples were collected and tested, which included 61 positive samples and 256 negative samples.

The Coronavirus Ag Rapid Test Cassette (Swab) results were compared to results of USFDA Emergency Use Authorized RT-PCR assays for SARS-CoV-2 in nasopharyngeal swab specimens. Overall study results are shown below.

Method		PCR		Total Results
Coronavirus Ag Rapid Test Cassette (Swab)	Results	Positive	Negative	Total Results
	Positive	59	2	61
	Negative	2	254	256
Total		61	256	317

#### \*Confidence Intervals

Relative Sensitivity: 96.72% (95%CI\*: 88.65%-99.60%)
Relative Specificity: 99.22% (95%CI\*: 97.21%-99.91%)
Accuracy: 98.74% (95%CI\*: 96.80%-99.66%)

## Benefits

- Rapid testing for SARS-CoV-2 antigen within the first ten days of symptom onset
- Rapid results within 15 minutes
- Facilitates patient treatment decisions quickly
- Simple, time-saving procedure
- No instrumentation required
- Very high sensitivity and specificity



# Ordering Information

Product DescriptionSpecimenCatalog No.FormatKit SizeCoronavirus Antigen Rapid TestNasopharyngeal SwabGCCOV-502a √Cassette20 Test/Kit

