



浙江东方基因生物制品股份有限公司  
Zhejiang Orient Gene Biotech Co.,LTD

## Manufacturer Statement

To whom it may concern,

The nucleocapsid phosphoprotein of SARS-CoV-2 was used as the mainly detected target in our Influenza&COVID-19 Ag Combo Rapid Test Cassette (Swab), COVID-19 Ag Rapid Test Cassette (Oral Fluid), Rapid COVID-19 Antigen Self-Test, Coronavirus Ag Rapid Test Cassette (Swab) these antigen tests can be used to detect both wild and mutant SARS-CoV-2 including the Alpha Variant (B.1.1.7), Beta Variant (B.1.351), Gamma Variant (P.1), Delta Variant (B.1.617.2), Lambda Variant (C.37), Omicron Variant (B.1.1.529).

However, these test reagents cannot distinguish these different strains.

Zhejiang orientgene Biotech Co.,Ltd.

2021.11.28



浙江省湖州市安吉县递铺街道阳光大道东段 3787 号

Address: 3787#, East Yangguang Avenue, Dipu Street, Anji 313300, Huzhou, Zhejiang, China.

电话 Tel: +86-572-5226111 传真 Fax: +86-572-5226222 邮编 P.C.: 313300